

INSTRUCTIONS FOR THE COLLECTION OF A URINE SPECIMEN

1. The preferred specimen for parasite analysis is the FIRST urine after arising from sleep for the day.
2. Have a Zip Lock bag or clean container handy that can hold the entire voided urine.

***** 3. Carefully open the specimen container. **Do not** pour out the liquid preservative.

4. Place an amount of urine into the transport container marked URINE to the bold red fill line on the label. Overfilling is not a problem.
5. Close the cap tightly.
6. Make sure your name and date of birth are on the label.
7. Place this transport container into the same bag with the stool sample and wrapped in an absorbent paper towel. Close the bag. Wash your hands and follow the shipping instructions.

Results

Upon completion of testing we will send you the results along with pertinent information and next steps by email and/or US mail. Thank you.